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PRECISION

ONE-Q-CEM in Dispenser

**Compomer cement
for permanent luting
in auto-mixing cartridges**

INSTRUCTIONS FOR USE

Description

ONE-Q-CEM in Dispenser is an auto-mixing, chemical cure cement based on compomer technology. It is radiopaque and easy to use.

ONE-Q-CEM in Dispenser flows easily and thus allows an unobjectionable positioning of restorations. Due to its translucency, highly esthetical results with ceramic materials are obtained.

As a compomer ONE-Q-CEM in Dispenser releases fluoride.

Composition

ONE-Q-CEM in Dispenser is made of a glass-ionomer inserted in a matrix of different dental resins, activator, catalyst, additives.

Indications

ONE-Q-CEM in Dispenser is indicated for the permanent luting of:

- All metal crowns and bridges including precious, semi-precious and non-precious metals
- Porcelain fused to metal crowns and bridges
- Gold inlays and onlays
- Porcelain fused to metal crowns with butt margins.

Advantages

- Correct mixing ratio of the components due to auto-mixing
- Fluoride release
- Radiopacity
- Flows easily to permit precise placement of the cement and easy placement of the restoration
- Exhibits a short setting time without high heat generation
- Insoluble in wet environment
- Insulator against thermal sensitivity
- Extremely low film thickness.

Side effects and contraindications

There are no known systemic side effects to date, however:

- Do not use ONE-Q-CEM in Dispenser in cases where known allergies to the components exist
- In some cases individual contact allergies may exist.

Precautions

Please use caution under all-ceramic crowns, inlays, onlays and veneers. The long term clinical stability with ceramics has not yet been fully established.

Avoid contact of ONE-Q-CEM in Dispenser under the secondary crowns when working with telescopic crowns and bridges.

Interactions with other substances

Do not use ONE-Q-CEM in Dispenser together with eugenol containing materials. We recommend the use of eugenol-free temporary cements.

Method of application

Preparation of the cartridge

Insert the cartridge into the specially developed dispenser for cements. Attach a new mixing tip for cements and turn it clockwise, until it locks into position. By pulling the trigger of the dispenser ONE-Q-CEM in Dispenser is pressed through the mixing tip and can be applied directly. Leave the mixing tip on the cartridge till the next use. It serves as a cap.

Preparation of the restoration

1. Sandblast (microetch) internal aspect of restoration with microetcher. High noble restorations do not need microetching. However, they can be tin-plated if desired.
2. Thoroughly wash and dry internal surface of restoration before cementation.
3. Place separating medium (e.g. Vaseline, glycerin) on external aspect of restoration for ease of clean up.

Preparation of the tooth

1. Cleanse prepared tooth surface thoroughly with pumice and water, or non-fluoride prophylaxis paste. Application of an antibacterial solution (i.e. chlorhexidine) is advantageous.
2. Rinse and dry with air syringe or cotton pellets, being careful to avoid pooling. Dentin should appear moist and shiny. Do not desiccate or contaminate.
3. To obtain the highest possible adhesion, use a primer and a bond.

Cementation of the restoration

1. Dispense ONE-Q-CEM in Dispenser either into the restoration or directly onto the restoration from the tip.

2. ONE-Q-CEM in Dispenser has a working time of 2 minutes.
3. Seat the restoration under light pressure. Verify complete seating, then have the patient bite slowly into occlusion. A bite stick may be used.
4. After placement, ONE-Q-CEM in Dispenser will set under conditions of the mouth within 1 minute and 30 seconds - 2 minutes.

Clean up

Remove all excess cement from around margins once cement has reached its soft, doughy state with a curette, explorer, or scaler. Begin with interproximal areas. Floss interproximal areas gently.

Important Notes

When using a cartridge for the first time please discard the first "pea" sized amount of extruded material.

Storage

Store ONE-Q-CEM in Dispenser at room temperature 18 to 23 °C (64 to 73 °F). Only use at room temperature (refrigerated material is more viscous and cures more slowly). Do not use after expiration date.

Packaging

1 double cartridge of Base and Catalyst of 100 g total, 10 mixing tips for cements.

For dental use only. Keep out of reach of children.

PRECISION products are carefully developed and extensively tested. We therefore guarantee their impeccable quality when used in accordance with their instructions for use. Application, use and processing of the product are beyond our control and are subject to the responsibility of the user.

For more information please contact:

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DENTALICA S.p.A. – Via Rimini 22 – I-20142 Milano – Italia