Simboli usati sulla confezione / Symbols used on packaging / Symboles utilisés sur les emballages / Símbolos utilizados en el embalaje / Auf der verpackung verwendete symbole / Symbole na opakowaniu / Simboluri utilizate pe ambalaj / Συμβόλα που χρησιμοποιούνται στη συσκευασία / Símbolos usados na embalagem

 Consultare le istruzioni d'uso / Consult instructions for use / Consulter le mode d'emploi / Consultar las instrucciones de uso / Die Gebrauchsanweisung lesen / Patrz: sposoby użycia / Consultați instrucțiunile de utilizare / Συμβουλευτείτε τις οδηγίες χρήσης / Consulte as instruções de uso	 Numero di catalogo / Catalogue number / Numéro de catalogue / Número de catálogo / Katalognummer / Numer katalogowy / Număr de catalog / Αριθμός καταλόγου / Número de catálogo
 Limiti di temperatura / Temperature limit / Limite de température / Limite de temperatura / Temperaturgrenze / Grænica temperatury / Limita de temperatură / Όριο θερμοκρασίας / Limite de temperatura	 Attention / Caution / Attention / Atención / Beachtung / Uwaga / Atentie / Προσοχή / Atenção
 Mantenere asciutto / Keep dry / Garder au sec / Manter seco / Bleib trocken / Utrzymuj suchotę / Păstrati uscat / Διατηρώ στεγνό / Manter seco	 Tenere lontano dalla luce solare / Keep away from sunlight / Tenir à l'écart de la lumière du soleil / Mantener alejado de la luz solar / Von direktem Sonnenlicht fernhalten / Chronić przed światłem słonecznym / A se păstă de lumina soarelui / Διατηρείται μακριά από την άμεση ηλιακή ακτινοβολία / Manter afstået da luz solar direta
 Data di scadenza / Use by date / Date limite d'utilisation / Utilizar según la fecha / Mindesthaltbarkeitsdatum / Termin przydatności do użycia / Data expirării / Τέλησι ημερομηνία ανάλυσης / Data de validade	 Rappresentante autorizzato / Authorized representative / Représentant autorisé / Representante autorizado / Bevollmächtigter Vertreter / Uprawniony przedstawiciel / Εξουσιοδοτημένος αντιπρόσωπος / Representante autorizado / Reprezentant autorizat
 Codice lotto / Batch code / Code de lot / Código de lote / Chargennummer / Kod partii / Cod lot / Καδώς αριθμός / Código de lote	 Produttore / Manufacturer / Fabricant / Fabricante / Hersteller / Producător / Κατασκευαστής / Fabricante

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NOSTRESS

EDTA lubricating cream for rotary files

INSTRUCTIONS FOR USE

NoStress is an EDTA cream with Urea Hydrogen Peroxide, used as an irrigation agent for chemo-mechanical preparation of the root canal.

Indications (use in conjunction with sodium hypochlorite)

- Debridement in the root canal.
- Dissolution of pulp remnants.
- Removal of smear layer.
- Facile movements of canal instruments.
- Microbicidal treatment within the canal.
- Clearance of dentin tubules and apical area.

Composition

- Ethylenediaminetetraacetic Acid (EDTA).
- Urea Hydrogen Peroxide (UHP).

Advantages and Effectiveness

- EDTA removes smear layer, softens dentin, and facilitates the removal of calcific obstruction.
- UHP reacts with sodium hypochlorite to produce oxygen bubbles that flush out debris within the canal.
- The lubricant medium facilitates instrument movement within the canal to avoid instrument binding or breakage.
- Unclogged dentin tubules and apical area permit better penetration of root canal medicaments and fillings.

Directions

1. Isolate the tooth to be treated with a rubber dam.
2. Open and clear the pulp chamber in a conventional manner.
3. Inject ample amount of NoStress into the pulp chamber or deposit NoStress directly on the canal shaping instrument.
4. Shape the root canal in a conventional manner.
5. Flush the root canal as needed with sodium hypochlorite solution.

6. Repeat step 3 ~ 5 until the canal is shaped and cleaned as desired.

Precautions

- For dental use only.
- Avoid extrusion of the excess cream out of the apex.
- Protective gloves are recommended.
- In case of accidental contact with eyes or prolonged exposure to oral tissues, flush immediately with ample amount of water.
- Do not use after expiry date.

Storage

- Store in a dry and cool place.
- Keep the cap tightly sealed after each use.
- Shelf life is 2 years.

Packaging

2 syringes of 6 g, 10 application tips.

For professional dental use only. Made in Korea.

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