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DETAIL FLUID NEXT

Addition curing silicone with intelligent matrix and high fluidity

INSTRUCTIONS FOR USE

Detail Fluid NEXt is a VPS precision impression material with low consistency for automatic mixing, based on addition-cured silicone.

Indication

Correction impression or double mixing technique in combina-tion with Detail Putty NEXtfor fabricating crowns, bridges, inlays and onlays.

Cartridge handling

The material is dispensed and mixed automatically by extruding through a mixing tip. When using a new cartridge remove the transport seal. Then attach

the mixing tip and, if required, an intraoral tip. Following use, the mixing tip should be used as a seal and should be left on the cartridge until the next applica-tion. The mixing tip should only be replaced with a new one when the cartridge is to be used again.

Recommended use

Preparing the Impression

Prepare the Detail Fluid NEXt wash material and check that the outlets of the safety cartridges are free from obstruction shortly before application. Used mixing tips serve as a seal until the next application. Never re-use the original seals of the safety cartridges. This can cause unintentional curing due to contami-nation of the base and catalyst paste. If retraction liquids have been used, thoroughly rinse the treatment site with water and dry it with oil-free air prior to taking the impression.

Correction Impression

Insert the impression tray loaded with Detail Putty NEXtinto the pa-tient's mouth. Maintain this position until the recommended time in mouth of the material. Thoroughly clean and dry the impres-sion after removal. Cut escape channels in the preliminary impression and remove high spots before taking the wash impression. Apply Detail Fluid NEXt wash material to the impression and syringe wash material around the prepared teeth, if neces-sary. Re-insert the loaded preliminary impression into the patient's mouth observing the working and setting times of the wash material. Do not wait longer than 30 minutes between the preliminary impression and wash impression to ensure that the two impression materials bond together.

Double mixing technique

While the assistant loads the tray, the operator can syringe around the prepared teeth with Detail Fluid NEXt. Load-ing the tray and syringing round the prepared teeth should be completed simultaneously if possible, as the wash material sets quicker than the material in the tray due to the higher intraoral temperature. Insert the impression tray into the patient's mouth. Maintain this position until the recommended time in mouth of the material.

Disinfecting and casting the impression

The impression can be disinfected with commercially available immersion solutions. We recommend using a 2% glutaralde-hyde solution. Ensure that the disinfectant is indicated for use with addition-cured impression silicones. Immerse the impression for approx. 10 to 15 minutes. Then rinse under running water for 15 seconds. Keep the impression at normal room temperature (max. 25°C / 77°F) and cast the impression no sooner than 30 minutes after removing it from the patient's mouth.*

The impression material is compatible with all commercial casts. We recommend using a special purpose hard plaster (ISO 6873, Type 3; e.g. Moldano, Bayer) for the cast or an ultra hard dental plaster for the forms (ISO 6873, Type 4:e.g. Fujirock II, GC). Impressions can be copper-plated or silver-plated.

** For further information see the comments of the DGZMK. Comments of the DGZMK 3/93 V 1.0; [www.dgzmk.de](#)*

ADA recommends following Council on Dental Materials, Instruments and Equipment. Disinfection of impressions. JADA 1991; 122(3):110; [www.ada.org](#)

Please note

- Do not use with condensation-cured or polyether impres-sion materials.
- Avoid contact of the impression with organic solvents or liquids containing solvent, as this can cause the material to swell and lead to inaccuracies.
- Retraction solutions should be thoroughly rinsed with water and dried with oil-free air before taking the impression to avoid them impairing the setting properties.
- Do not use with latex gloves.
- In rare cases it cannot be ruled out that predisposed pa-tients may develop hypersensitivity to individual compo-nents of the material.

Technical data Detail Fluid NEXt

Working time (incl. mixing time) ≤ 1:15 min.

Recommended time in mouth ≥ 2:00 min.

Strain in compression approx. ~ 3,8%

Compression set ≤ 0,4%

Linear change in length (24 h) ≤ 0,35%

Times refer to normal room temperature (23°C / 73°F) and a normal relative humidity (50%). Higher temperatures shorten, lower temperatures increase those times.

Composition

Vinyl polysiloxanes (VPS), silicon dioxide, platinum catalyst, ten-sides and in small amounts stabilizers and additives.

Classification

DIN EN ISO 4823:2000, Type 3, low consistency- light bodied.

Storage

Store at room temperature (15-25°C / 59-77°F). Do not use after the expiration date.

Available product sizes

Detail Fluid NEXt: 2 cartridges @ 50 ml, 12 Mixing tips, 12 Intraoral tips

For professional dental use only. Keep away from children.

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