NOVABIO

(POLIVINILIDENFLUORURO)

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Synthetic sterile non-absorbable sutures

DESCRIPTION

NOVABIO is a monofilament, synthetic, sterile, non-absorbable suture for surgical use, made of polyvinylidene fluoride. NOVABIO is colored with Phthalocyanine blue. NOVABIO sutures are available in a range of gauge sizes and lengths, non-needed or attached to stainless steel needles of different types and sizes. The needles' surface is provided with a silicone coat to make the tissues' penetration easier and less traumatic for the patient. NOVABIO sutures are available with PTFE pledgets. NOVABIO complies with the requirements of the European Pharmacopoeia's monograph "Sterile, non-absorbable sutures". NOVABIO is natural rubber latex free.

Indications

NOVABIO is recommended for its long resistance over time and its knots' stability, especially for long term implantation in vascular surgery. It is usually employed for sutures and ligatures in soft tissues, including use in neurosurgery, ophthalmic and cardiovascular surgery. Moreover, NOVABIO is used on the skin, in particular in plastic surgery.

Applications

The suture has to be chosen depending on the patient's conditions, the surgeon's experience, the surgical technique and the wound's size.

Properties

NOVABIO elicits a minimal initial inflammatory reaction followed by a slow and progressive encapsulation of connective tissue.

Contraindications

Not Known.

Warnings/Precautions/Interactions

NOVABIO sutures must be exclusively used by specialized personnel, qualified for the surgical practice. The correct application procedure minimizes the risk of damages. The surgeon, who is trained for selecting the product and the most proper format to minimize the risk of damages, has the responsibility to choose the correct format to employ depending on the situation. NOVABIO's available commercial variants satisfy different needs.

Before using NOVABIO, the user must verify the integrity of the package.

Open the primary package just before the use, respecting asepsis' rules. While removing the package, then, the user must be sure to preserve the sterility of the product contained in it. Take out the inner sterile package wearing sterile gloves.

The user must have a good knowledge of surgical techniques in which NOVABIO is employed.

As any suture's material, pay attention not to damage the wire while handling it. In particular, avoid bending or crushing with surgical instruments such as needle-holders. While tightening the suture, pull the wire between the needle and the needle's penetration channel. Do not pull the wire with too much strength or over sharp objects. While pulling the wire, avoid friction with gloves to prevent the abrasion of the surface and damages on the structure of the wire itself.

To avoid damages, the needle has to be grabbed between 1/3 and 1/2 of the distance from the attachment end to the point. Avoid to deform the needle because it will induce a loss in stability of the needle itself. While handling surgical needles, the user must pay attention not to get hurt, being them a possible source of infection. Use needles must be disposed following the related directive (to avoid a possible risk of infection).

The risk of wound's dehiscence depends on the localization of the wound and the suture employed. For this reason, the user must have a good knowledge of surgical techniques in which NOVABIO is employed. NOVABIO's prolonged contact with salt solutions may induce, as any foreign body, the formation of concretions (urinary and biliary tracts).

NOVABIO, when it is recommended, is employed in usual suturing and knotting surgical techniques depending on the personal experience of the user. Used needles must be treated following the safety procedures for the disposal of hospital wastes.

Side effects

Side effects associated with this product are represented by an initial tissue inflammatory reaction of small entity and a temporary local irritation of the wound. As any foreign body, NOVABIO suture may aggravate a pre-existent infection.

Sterility

NOVABIO sutures are sterilized with Ethylene Oxide and must not be re-sterilized. Open or damaged packages must not be used.

Storage

It is recommended to store the device in a dry place, away from light sources. Do not use after the expiration date.

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