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BIOSILK
(VIRGIN SILK)
STERILE NATURAL NON-ABSORBABLE SUTURE**DESCRIPTION**

BIOSILK is a sterile, non-absorbable suture for surgical use composed by an organic protein called fibroin; this protein is obtained from silkworms' larvae (*Bombyx Mori*) of the *Bombycidae* family. BIOSILK is available in a range of gauge sizes and lengths, non-needled or attached to stainless steel needles of different types and sizes.

Needles' surface is provided with a silicone coat to make tissue penetration easier and less traumatic for the patient. BIOSILK complies with the European Pharmacopoeia's monograph for Sterile Braided Silk sutures. Each fiber, used to manufacture braided wires, is previously treated to remove natural waxes and gum resins. Successively, after braiding, wires are coated with a special

composition of waxes and silicone elastomers. Braided silk is available in black (C.I. 53185). As regards Virgin Silk, the sericin resin is not removed in order to keep together a certain number of braided silk wires. The Virgin Silk is available in Blue, colored with methylene blue dye (C.I. 52015). BIOSILK is natural rubber latex free.

INDICATIONS

BIOSILK is recommended for the approximation and/or ligature of soft tissues, including cardiovascular, ophthalmic, neurosurgical and oral surgery.

APPLICATIONS

Sutures should be selected and implanted depending on patient conditions, surgical experience, surgical technique and wound size.

PROPERTIES

BIOSILK elicits an inflammation reaction in the tissue, after which it is gradually encapsulated by fibrous connective tissue. Although BIOSILK is not absorbed, an in vivo biodegradation of the proteinaceous silk fibres can lead to long-term gradual loss of its tensile strength.

CONTRAINDICATIONS

This suture is not recommended on patients sensitive or allergic to silk.

WARNINGS/ PRECAUTIONS/ INTERACTIONS

Before using BIOSILK to close the wound, the user must have developed practice and experience with surgeries and surgical techniques which employ non-absorbable sutures: the risk of wound’s dehiscence, in fact, depends on the localization of the wound and the suture employed. As any foreign body, prolonged contact of the suture with salt solutions, such as those of urinary and biliary tracts, may induce the formation of calculi. While treating infected or contaminated wounds, it is necessary to employ the proper surgical practice. Before using BIOSILK, the user must verify the integrity of the package. Open the primary package just before the use, respecting asepsis’ rules. While removing the package, then, the user must be sure to preserve the sterility of the product contained in it. Take out the inner sterile package wearing sterile gloves. While using BIOSILK, as any other suture, pay attention not to damage the suture itself, avoiding tightening and/or crushing wires with forceps, needle-holders and other surgical instruments. As any suture, in order to guarantee the safety of the knot, it is recommended to apply the surgical technique of the flat and squared knot, with eventual additional knots, depending on the surgical circumstance and the user’s experience. While handling surgical needles, pay attention not to damage them. The needle must be clamped in the needle-holder between one third and one half of the distance from the attachment end to the point. If it is grabbed too near to the tip, it will obstruct the penetration and/or induce its breaking. If it is grabbed too near to the needle’s eye, instead, it might bend or break. Trying to modify the curvature of the needle may reduce its bending and breaking resistance. While handling surgical needles, pay attention not to get hurt. Used needles must be disposed in proper containers.

SIDE EFFECTS

Side effects associated with this product include an allergic reaction in patients with a known sensitivity to silk and an initial tissue inflammatory reaction or a temporary irritation of the wound. As any foreign body, silk sutures might aggravate a pre-existent infection of the wound.

STERILITY

BIOSILK sutures are sterilized with Ethylene Oxide and must not be re-sterilized. Do not use if the package is open or damaged.

STORAGE

It is recommended to store the device in a dry place, away from light sources. Do not use after the expiration date.

LEGEND OF THE SYMBOLS ON THE PACKAGE

	The product has to be stored in a dry place
	Keep the product away from the light
	Do not use if the package is damaged
	Carefully read the instructions for use
	Do not re-use
	Do not re-sterilize
	Reference code
	Lot number
	Expiration date (use within the date indicated, month and year)
	Manufacturer (Bio Sud Medical Systems S.r.l.)
	Sterile product. Sterilization method: Ethylene Oxide
	Colored, multifilament, non-absorbable suture
	Needle’s typology: curvature, dimension and shape of the tip (an example is shown)
	Suture size following the United States Pharmacopoeia (USP) and the European Pharmacopoeia (EP – Gauge number)
	Suture’s length
	QR code, reporting the reference to the product code, the lot number and the expiry date of the medical device
	CE marking and identification number of the notified organism. The product complies with the essential requirements of the Medical Devices Directive 93/42/CEE

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Producător (Bio Sud Medical Systems S.r.l.)