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**POLYBIO
(POLYESTER)
SYNTHETIC STERILE NON-ABSORBABLE SUTURE**

DESCRIPTION

Polybio is a braided, sterile, non- absorbable suture for surgical use made of coated poly-ethylene terephthalate.

The wire is coated with a silicone elastomer MED2174, that acts as lubricant and make the passage through tissues mechanically easier and fluent.

The wire is available in white and green, colored with D&C Green N° 6 dye (C. I. 61565).

POLYBIO is available in a range of gauge sizes and lengths, non-neededled or attached to stainless steel needles of different types and sizes. Needles' surface is provided with a silicone coating to make the tissues' penetration easier and less traumatic for the patient. POLYBIO sutures are available with PTFE pledgets. The sutures are natural rubber latex free.

Polybio complies with the requirements of the European Pharamacopoeia's monograph "Sterile, non-absorbable sutures".

INDICATIONS

Polybio is indicated for suture and ligature in soft tissues, for cardiovascular, ophthalmic surgeries and neurosurgery.

APPLICATIONS

Sutures should be selected and implanted depending on patient conditions, surgical experience, surgical technique and wound size.

PROPERTIES

Polybio sutures elicit a minimal initial inflammatory reaction in tissues and ingrowth of fibrous connective tissue. Generally, polyester doesn't show any relevant loss in resistance over time.

CONTRAINDICATIONS

Not known.

WARNINGS/ PRECAUTIONS/ INTERACTIONS

Polybio sutures must be exclusively used by specialized personnel, qualified for the surgical practice. The application of the correct procedure minimizes the risk of damages. The surgeon, who is trained to select the proper product and format in order to minimize the risk of damages, has the responsibility to choose the correct format to employ in different situations. Polybio available commercial variants satisfy different needs. Before using Polybio, the user must verify the integrity of the package. Open the primary package just before the use, respecting asepsis' rules. While removing the package, then, the operator must be sure to preserve the sterility of the product contained in it.

Take out the inner sterile package while wearing sterile gloves.

Before using Polybio for suturing wounds, the operator must have experience with surgical procedures and with non-absorbable sutures. Dehiscence's risk, in fact, depends on the localization of the wound and the suture employed.

when it is recommended, it is used for the usual surgical techniques of suturing and knotting according to the user's personal experience.

As any foreign body, prolonged contact of a suture with salt solutions, such as those of the urinary and biliary tracts, may induce the formation of calculi. It is necessary to follow the proper surgical practice for treating infected or contaminated wounds.

As any suture, pay attention not to damage the wire while handling it. In particular, avoid to crimp or crush the wire with surgical instruments such as needle-holders. While tightening the suture, pull the wire between the needle and the needle's penetration channel. Do not pull the wire with too much strength or over sharp objects. Avoid the friction with gloves in order to prevent the risk of superficial abrasion and damages to the wire's structure itself.

In order not to damage the needle, it has to be grabbed between 1/3 and 1/2 of the distance from the attachment end to the point. Avoid to deform the needle because it will lead to a loss in stability of the needle itself. While using surgical needles, the user must pay attention not to get hurt because they are a possible source of infection. Used needles must be treated respecting the safety procedures for the disposal of hospital wastes.

SIDE EFFECTS

Side effects associated with this product are represented by an initial tissue inflammatory reaction of small entity and a temporary local irritation of the wound. As any foreign body, the suture may aggravate a pre-existent infection.

STERILITY

Polybio sutures are sterilized with Ethylene Oxide and must not be re-sterilized. Do not use if the package is open or damaged.

STORAGE

It is recommended to store the device in a dry place, away from light sources. Do not use after the expiration date.

LEGEND OF THE SYMBOLS ON THE PACKAGE

	Do not use if the package is damaged
	Carefully read the instructions for use
	Do not re-use
	Do not re-sterilize
	Reference code
	Lot number
	Expiration date (use within the date indicated, month and year)
	Manufacturer (Bio Sud Medical Systems S.r.l.)
	Sterile Product. Sterilization method: Ethylene Oxide
	Keep dry
	Keep the product away from the light
	Multifilament, braided, undyed, coated, non-absorbable suture
	Multifilament, braided, colored, coated, non-absorbable suture
	Pre-assembled suture with pledget
	Needle's typology: curvature, size and shape of the tip (an example is shown)
USP __ EP __	Suture size following United States Pharmacopoeia (USP) and European Pharmacopoeia (EP – Gauge number)
__ cm	Suture's length
	Data-matrix, reporting the reference to the product code, the lot number and the expiry date of the medical device
	CE marking and identification number of the notified Organism. The product complies with the essential requirements of the Medical Devices Directive 93/42/CEE

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Producător (Bio Sud Medical Systems S.r.l.)