Manipoli per scaling

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Scaling Handpieces

Xpedent Handpiece series, which has 4 specifications, is mainly used for ultrasonic scaling devices. Please confirm whether the purchased parts can match the equipment.

Handpieces must be sterilised before first use and after each use (they are not delivered sterile).

Cleaning:

The removal of organic deposits is very important and will achieve best results by putting them immediately after use into 70% Isopropyl alcohol (Isopropanol) or a low alcohol cleaning solution.

After that use an ultrasonic bath with a low foaming detergent. Alternatively, they can be cleaned in warm water with a low foaming detergent and an autoclavable brush.

After the cleaning:

Rinse the cleaned handpieces thoroughly in warm water to remove all chemicals and let them dry completely. Therefore, the use of an antirust agent before steam sterilisation is acceptable, although it is not required. Pack them in suitable packaging, bags, trays or cassettes.

Check them with spore test sets. Biological monitoring devices should be used at least once a week. In addition, every time if new packaging, new staff, new procedures, new equipment or repaired equipment is used, a biological monitoring should be added to check if these changes are still effective.

Warning:

Do not use closed containers during sterilisation with steam or ethylene oxide. The cleaning and Disinfecting device must be approved by its manufacturer for cleaning and disinfection be released by dental instruments and EN ISO 15883 1

Sterilise with ethylene oxide

Put the dried packaged materials in the unit and work according to the manufacturer's instructions.

Sterilise with steam

Place the dried wrapped tip in an autoclave. Recommended settings are 134 degrees for 18 minutes.

Regarding sterilisation:

Ethylene oxide is the preferred method for sterilisation. Steam autoclaving is the preferred method of dentistry. Dry heat should NOT be used for materials which can be sterilised by steam.

Chemical agents are not recommended since their effectiveness cannot be monitored biologically, their cycle times are longer (6-10 hours), and all rinsing and handling must be aseptic.

Dry Heat and chemical steam sterilisation are considered ineffective methods, with the risk of damaging materials.

- The minimum sterilisation cycle for handpiece and torque wrench is 250 cycles.

- The handpiece must not be soaked in solution. It should be wiped clean with a gauze or soft sponge using 70% isopropanol and detergent or appropriate disinfection equivalent. Autoclaving is recommended to avoid cross patient infection.
- The cover can be wiped off with a soft cloth. Use 70% isopropanol and detergent or appropriate disinfection equivalent.

Signs	Notes on the signs
	CE mark
	134°C Maximum autoclave temperature
	Manufacturer
	Medical Device
	Serial Number
	European Representative

Selected, distributed and guaranteed by

DENTALICA S.p.A. Via Rimini 22 | 20142 Milano - Italia
Tel. +39 02 89.598.1 | Fax +39 02 89.50.42.49
www.dentalica.com | dentalica@dentalica.com

Guilin Yikeshi Medical Instrument Co., Ltd
D08 High Tech Park Chaoyang Road
Guilin, 541004 PRC

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Xpedent Austria GmbH